



SHAMIR

SHAMIR QUICK REFERENCE GUIDE

			Hard Resin	Transitions®	Transitions® XTRActive™	Polarized	Drivewear®	DLC™ (Trivex)	DLC™ (Trivex) Transitions®	DLC™ (Trivex) Transitions® XTRActive™	DLC™ (Trivex) Polarized	1.56	Polycarbonate	Polycarbonate Transitions®	Polycarbonate Transitions® XTRActive™	Polycarbonate Polarized	Polycarbonate Drivewear®	Superlite™ 1.6	Superlite 1.6 Transitions®	Superlite 1.6 Polarized	Superlite 1.67	Superlite 1.67 Transitions® XTRActive™	Superlite 1.67 Polarized	Superlite 1.74				
	Fitting Height	INDEX	1.498	1.5	1.5	1.498	1.5	1.53	1.53	1.53	1.53	1.56	1.586	1.586	1.586	1.586	1.586	1.6	1.6	1.6	1.67	1.67	1.67	1.67				
Shamir Autograph FreeFrame™		Power Range	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-14.00 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	-14.75 to +8.25	-14.75 to +8.25	-14.75 to +8.25	-14.75 to +8.25	-18.75 to +15.50			
		Nominal Base Curve*	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7, 8	0.5, 2, 3, 4, 5, 6, 7, 8	0.5, 2, 3, 4, 5, 6, 7, 8	0.5, 2, 3, 4, 5, 6, 7, 8	0.5, 1.5, 2.5, 3.5, 4.75, 6, 7.25, 8.25, 9, 9.75, 10.5, 11.25, 12			
		Add	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00	0.75 to 4.00		
		Color	Clear	Gray & Brown	Gray	Gray & Brown	Gray	Drivewear®	Clear	Gray & Brown	Gray	Gray & Brown	Clear	Clear	Gray & Brown	Gray	Gray & Brown	Drivewear®	Clear	Gray & Brown	Gray & Brown	Clear	Gray & Brown	Gray & Brown	Gray & Brown	Clear		
Shamir Autograph Attitude™		Power Range	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00	15mm: -6.50 to +5.25 18mm: -6.50 to +6.00		
		Nominal Base Curve*	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8, 9, 10	15mm: 6, 7, 8, 9, 10 18mm: 6, 7, 8, 9, 10	15mm: 6, 7, 8, 9, 10 18mm: 6, 7, 8, 9, 10	15mm: 6, 7, 8, 9, 10 18mm: 6, 7, 8, 9, 10	15mm: 6, 7, 8, 9, 10 18mm: 6, 7, 8, 9, 10	15mm: 6, 7, 8, 9, 10 18mm: 6, 7, 8, 9, 10	15mm: 6, 7, 8, 9, 10 18mm: 6, 7, 8, 9, 10	15mm: 6, 7, 8, 9, 10 18mm: 6, 7, 8, 9, 10	15mm: 6, 7, 8, 9, 10 18mm: 6, 7, 8, 9, 10	15mm: 6, 7, 8, 9, 10 18mm: 6, 7, 8, 9, 10	15mm: 6, 7, 8, 9, 10 18mm: 6, 7, 8, 9, 10	15mm: 6, 7, 8, 9, 10 18mm: 6, 7, 8, 9, 10	15mm: 6, 7, 8, 9, 10 18mm: 6, 7, 8, 9, 10	15mm: 6, 7, 8, 9, 10 18mm: 6, 7, 8, 9, 10			
		Add	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50
		Color	Clear	Gray & Brown	Gray	Gray & Brown	Drivewear®	Clear	Gray & Brown	Gray	Gray & Brown	Gray	Clear	Clear	Gray & Brown	Gray	Gray & Brown	Drivewear®	Clear	Gray & Brown	Gray & Brown	Clear	Gray & Brown	Gray & Brown	Gray & Brown	Clear		
Shamir Autograph Plus™		Power Range	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-14.00 to +6.50	-14.00 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	-14.75 to +8.00	-14.75 to +8.00	-14.75 to +8.00	-14.75 to +8.00	-18.75 to +14.75		
		Nominal Base Curve*	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 1.5, 2.5, 3.5, 4.75, 6, 7.25, 8.25, 9, 9.75, 10.5, 11.25, 12		
		Add	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	
		Color	Clear	Gray & Brown	Gray	Gray & Brown	Drivewear®	Clear	Gray & Brown	Gray	Gray & Brown	Gray	Clear	Clear	Gray & Brown	Gray	Gray & Brown	Drivewear®	Clear	Gray & Brown	Gray & Brown	Clear	Gray & Brown	Gray & Brown	Gray & Brown	Clear		
Shamir Spectrum™		Power Range	-7.00 to +4.00	-7.00 to +4.00	-7.00 to +4.00	-7.00 to +4.00	-7.00 to +4.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-12.00 to +8.00	-12.00 to +8.00	-12.00 to +8.00	-12.00 to +8.00	-12.00 to +8.00		
		Nominal Base Curve*	2, 3, 7.5, 4.75, 5.25, 6, 7, 8	2, 3, 7.5, 4.75, 5.25, 6, 7, 8	2, 3, 7.5, 4.75, 5.25, 6, 7, 8	2, 3, 7.5, 4.75, 5.25, 6, 7, 8	2, 3, 7.5, 4.75, 5.25, 6, 7, 8	2, 3, 7.5, 4.75, 5.25, 6, 7, 8	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	2, 3, 4, 5, 6, 7, 8, 9	2, 3, 4, 5, 25, 6, 25, 7.5, 8.5, 9	2, 3, 4, 5, 25, 6, 25, 7.5, 8.5, 9	2, 3, 4, 5, 25, 6, 25, 7.5, 8.5, 9	2, 3, 4, 5, 25, 6, 25, 7.5, 8.5, 9	2, 3, 4, 5, 25, 6, 25, 7.5, 8.5, 9	2, 3, 4, 5, 6, 7, 8, 9	2, 3, 4, 5, 6, 7, 8, 9	2, 3, 4, 5, 6, 7, 8, 9	2, 3, 4, 5, 6, 7, 8, 9	2, 3, 4, 5, 6, 7, 8, 9	2, 3, 4, 5, 6, 7, 8, 9	2, 3, 4, 5, 6, 7, 8, 9	2, 3, 4, 5, 6, 7, 8, 9		
		Add	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	0.75 to 3.50	
		Color	Clear	Gray & Brown	Gray	Gray & Brown	Drivewear®	Clear	Gray & Brown	Gray	Gray & Brown	Gray	Clear	Clear	Gray & Brown	Gray	Gray & Brown	Drivewear®	Clear	Gray & Brown	Gray & Brown	Clear	Gray & Brown	Gray & Brown	Gray & Brown	Clear		
Shamir Computer™ & Workspace™		Power Range	-6.00 to +6.00										-6.00 to +8.00						-6.00 to +8.00									
		Nominal Base Curve*	2, 3, 7.5, 5.25, 6.50, 8											2, 3, 5, 6.25, 8, 9.5					2, 3, 5, 6.25, 8, 9.5									
		Add	0.75 to 3.50											0.75 to 3.50					0.75 to 3.50									
		Color	Clear											Clear						Clear								
Computer 14mm	Invisible Marking	15	FFC	14									PC	FFC	14													
		15	FFW	14										PC	FFW	14												

\* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



SHAMIR

SHAMIR QUICK REFERENCE GUIDE

			Hard Resin	Transitions®	Transitions® XTRActive™	Polarized	Drivewear®	DLC™ (Trivex)	DLC™ (Trivex) Transitions®	DLC™ (Trivex) Transitions® XTRActive™	DLC™ (Trivex) Polarized	1.56	Polycarbonate	Polycarbonate Transitions®	Polycarbonate Transitions® XTRActive™	Polycarbonate Polarized	Polycarbonate Drivewear®	Superlite™ 1.6	Superlite 1.6 Transitions®	Superlite 1.6 Polarized	Superlite 1.67	Superlite 1.67 Transitions®	Superlite 1.67 Transitions® XTRActive™	Superlite 1.67 Polarized	Superlite 1.74	
	Fitting Height	INDEX	1.498	1.5	1.5	1.498	1.5	1.53	1.53	1.53	1.53	1.56	1.586	1.586	1.586	1.586	1.586	1.6	1.6	1.6	1.67	1.67	1.67	1.67	1.74	
Shamir Relax™	Power Range		-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00								-14.00 to +8.00					-14.00 to +8.00								
	Nominal Base Curve*		0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8								0.5, 2, 3, 4, 2.5, 5.25, 6.25, 7.5, 8.5					0.5, 2, 3, 4, 5, 6, 7, 8								
	Add		0.75 to 3.50	0.75 to 3.50	0.75 to 3.50								0.75 to 3.50					0.75 to 3.50								
	Color		Clear	Gray & Brown	Gray								Clear					Clear								
Shamir FirstPAL™	Power Range		-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00							-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00		-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00						
	Nominal Base Curve*		0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8							0.5, 2, 3, 4, 5, 2.5, 6.25, 7.5	0.5, 2, 3, 4, 5, 2.5, 6.25, 7.5	0.5, 2, 3, 4, 5, 2.5, 6.25, 7.5	0.5, 2, 3, 4, 5, 2.5, 6.25, 7.5		0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7	0.5, 2, 3, 4, 5, 6, 7						
	Add		0.75 to 1.50	0.75 to 1.50	0.75 to 1.50	0.75 to 1.50							0.75 to 1.50	0.75 to 1.50	0.75 to 1.50	0.75 to 1.50		0.75 to 1.50	0.75 to 1.50	0.75 to 1.50						
	Color		Clear	Gray & Brown	Gray	Gray & Brown							Clear	Gray & Brown	Gray	Gray & Brown		Clear	Gray & Brown	Gray & Brown						
Shamir Cool™	Power Range		-12.25 to +6.00					-13.25 to +6.00					-14.75 to +6.00					-15.00 to +6.00								
	Nominal Base Curve*		0.5, 2, 3, 7.5, 5.25, 6, 7, 8					0.5, 2, 3, 4, 5, 6, 7					0.5, 2, 3, 4, 5.25, 6.25, 7.5					0.5, 2, 3, 4, 5, 6, 7								
	Add		0.75 to 3.50					0.75 to 3.50					0.75 to 3.50					0.75 to 3.50								
	Color		Clear					Clear					Clear					Clear								
Shamir Cool SV™	Power Range		-12.25 to +6.00					-14.00 to +7.50					-14.00 to +8.75					-14.00 to +8.00								
	Nominal Base Curve*		0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8, 10					0.5, 2, 3, 4, 5, 6, 7, 8, 9					0.5, 2, 3, 4, 5.25, 6.25, 7.5, 8.5					0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75								
	Add		0.75 to 3.50					0.75 to 3.50					0.75 to 3.50					0.75 to 3.50								
	Color		Clear					Clear					Clear					Clear								
Shamir Attitude Cool SV™	Power Range		-6.00 to +6.00					-6.00 to +6.00					-6.50 to +6.00					-6.50 to +6.00								
	Nominal Base Curve*		6, 7, 8					6, 7, 8					4, 6, 7, 8					4, 6, 7, 8, 9, 10								
	Add		Clear & Polarized					Clear & Polarized					Clear & Polarized					Clear & Polarized								
	Color		Clear & Polarized					Clear & Polarized					Clear & Polarized					Clear & Polarized								
Attitude Cool SV™	Power Range		-6.00 to +6.00					-6.00 to +6.00					-6.50 to +6.00					-6.50 to +6.00								
	Nominal Base Curve*		6, 7, 8					6, 7, 8					4, 6, 7, 8					4, 6, 7, 8, 9, 10								

\* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



SHAMIR

SHAMIR QUICK REFERENCE GUIDE

			Hard Resin	Transitions®	Transitions® XTRActive™	Polarized	Drivewear®	DLC™ (TriVex)	DLC™ (TriVex) Transitions®	DLC™ (TriVex) Transitions® XTRActive™	DLC™ (TriVex) Polarized	1.56	Polycarbonate	Polycarbonate Transitions®	Polycarbonate Transitions® XTRActive™	Polycarbonate Polarized	Polycarbonate Drivewear®	Superlite™ 1.6	Superlite 1.6 Transitions®	Superlite 1.6 Polarized	Superlite 1.67	Superlite 1.67 Transitions®	Superlite 1.67 XTRActive™	Superlite 1.67 Polarized	Superlite 1.74	
	Fitting Height	INDEX	1.498	1.5	1.5	1.498	1.5	1.53	1.53	1.53	1.53	1.56	1.586	1.586	1.586	1.586	1.586	1.6	1.6	1.6	1.67	1.67	1.67	1.67	1.74	
Shamir Ski™	Power Range							-5.00 to +4.00					15mm: -6.50 to +6.00 18mm: -6.50 to +5.25					15mm: -6.50 to +6.00 18mm: -6.50 to +5.25							-6.00 to +5.00	
	Nominal Base Curve*							6, 7, 8					15mm: 4, 5, 6, 7, 8, 9, 10 18mm: 4, 5, 6, 7, 8					15mm: 4, 5, 6, 7, 8, 9, 10 18mm: 4, 5, 6, 7, 8						15mm: 6, 7, 8 18mm: 5, 6, 7, 8		
	Add							0.75 to 3.50					0.75 to 3.50					0.75 to 3.50								0.75 to 3.50
	Color							Clear					Clear					Clear				Clear				
Shamir Ski SV™	15mm	Invisible Marking 15mm						TX ski5					PC ski5					16 ski5				67 ski5				
	18mm	Invisible Marking 18mm						TX ski8					PC ski8					16 ski8				67 ski8				
	Power Range							-6.00 to +6.00					-6.50 to +6.00					-6.50 to +6.00								-6.00 to +6.00
	Nominal Base Curve*							6, 7, 8					4, 6, 7, 8					4, 6, 7, 8, 9, 10								6, 7, 8
Shamir Golf™	Power Range							-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25
	Nominal Base Curve*							6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8
	Add							1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	1.00 to 3.50	
	Color							Clear	Gray & Brown	Gray	Gray & Brown	Clear	Gray & Brown	Gray	Gray & Brown	Gray & Brown	Clear	Gray & Brown	Gray & Brown	Gray & Brown	Clear	Gray & Brown	Gray & Brown	Gray & Brown	Gray & Brown	Gray & Brown
Shamir Golf SV™	19mm	Invisible Marking 19mm						TX SV	TTX SV	XX SV	PTX SV	56 SV	PC SV	PT SV	PX SV	PPZ SV	PW SV	16 SV	T6 SV	PZ6 SV	67 SV	T67 SV	X67 SV	PZ67 SV	74 SV	
	Power Range							-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	
	Nominal Base Curve*							6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8
	Color							Clear	Gray & Brown	Gray	Gray & Brown	Clear	Gray & Brown	Gray	Gray & Brown	Gray & Brown	Clear	Gray & Brown	Gray & Brown	Gray & Brown	Clear	Gray & Brown	Gray & Brown	Gray & Brown	Gray & Brown	
Shamir Smart Attitude™	Power Range		-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00
	Nominal Base Curve*		6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	4, 6, 7, 8	4, 6, 7, 8	4, 6, 7, 8	4, 6, 7, 8	4, 6, 7, 8	4, 6, 7, 8, 9, 10	4, 6, 7, 8, 9, 10	4, 6, 7, 8, 9, 10	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8	6, 7, 8
	Color		Clear	Gray & Brown	Gray	Gray & Brown	Drivewear®	Clear	Gray & Brown	Gray	Gray & Brown	Clear	Clear	Gray & Brown	Gray	Gray & Brown	Drivewear®	Clear	Gray & Brown	Gray & Brown	Clear	Gray & Brown	Gray & Brown	Gray & Brown	Gray & Brown	Gray & Brown
Shamir Smart SV™	Invisible Marking		15 FF SA	TR FF SA	X5 FF SA	PZ FF SA	DW FF SA	TX FF SA	TTX FF SA	XX FF SA	PTX FF SA	56 FF SA	PC FF SA	PT FF SA	PX FF SA	PPZ FF SA	PW FF SA	16 FF SA	T6 FF SA	PZ6 FF SA	67 FF SA	T67 FF SA	X67 FF SA	PZ67 FF SA	74 FF SA	
	Power Range		-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +7.50	-14.00 to +7.50	-14.00 to +7.50	-14.00 to +7.50	-14.00 to +7.50	-14.00 to +8.50	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-16.00 to +16.00	-16.00 to +16.00	-16.00 to +16.00	-16.00 to +16.00	-16.00 to +16.00	-18.75 to +15.50
	Nominal Base Curve*		0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8, 10	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8, 10	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8, 10	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8, 10	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8, 10	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 25, 6.25, 7.5, 8.5	0.5, 2, 3, 4, 5, 25, 6.25, 7.5, 8.5	0.5, 2, 3, 4, 5, 25, 6.25, 7.5, 8.5	0.5, 2, 3, 4, 5, 25, 6.25, 7.5, 8.5	0.5, 2, 3, 4, 5, 25, 6.25, 7.5, 8.5	0.5, 2, 3, 4, 5, 25, 6.25, 7.5, 8.5	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75, 9.5, 11, 12.5	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75, 9.5, 11, 12.5	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75, 9.5, 11, 12.5	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75, 9.5, 11, 12.5	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75, 9.5, 11, 12.5

\* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



SHAMIR

SHAMIR QUICK REFERENCE GUIDE

			Hard Resin	Transitions®	Transitions® XTRActive™	Polarized	Drivewear®	DLC™ [Triplex]	DLC™ [Triplex] Transitions®	DLC™ [Triplex] Transitions® XTRActive™	DLC™ [Triplex] Polarized	1.56	Polycarbonate	Polycarbonate Transitions®	Polycarbonate Transitions® XTRActive™	Polycarbonate Polarized	Polycarbonate Drivewear®	Superlite™ 1.6	Superlite 1.6 Transitions®	Superlite 1.6 Polarized	Superlite 1.67	Superlite 1.67 Transitions®	Superlite 1.67 XTRActive™	Superlite 1.67 Polarized	Superlite 1.74
	Fitting Height	INDEX	1.498	1.5	1.5	1.498	1.5	1.53	1.53	1.53	1.53	1.56	1.586	1.586	1.586	1.586	1.586	1.6	1.6	1.6	1.67	1.67	1.67	1.67	1.74
Shamir Smart SV As Worn™	Power Range		-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +7.50	-14.00 to +7.50	-14.00 to +7.50	-14.00 to +7.50	-14.00 to +8.50	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-16.00 to +16.00	-16.00 to +16.00	-16.00 to +16.00	-16.00 to +16.00	-20.00 to +15.50
	Nominal Base Curve*		0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8, 10	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8, 10	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8, 10	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8, 10	0.5, 2, 3, 7.5, 4.25, 5.25, 6, 7, 8, 10	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75, 9.5, 11, 12.5	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75, 9.5, 11, 12.5	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75, 9.5, 11, 12.5	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75, 9.5, 11, 12.5	0.5, 1, 1.5, 2.5, 3.5, 4.75, 6, 7.25, 8.25, 9, 9.75, 10.5, 11.25, 12
	Color	Clear	Gray & Brown	Gray	Gray & Brown	Drivewear®	Clear	Gray & Brown	Gray	Gray & Brown	Clear	Clear	Gray & Brown	Gray	Gray & Brown	Drivewear®	Clear	Gray & Brown	Gray & Brown	Clear	Gray & Brown	Gray	Gray & Brown	Clear	Clear
	Invisible Marking	15 FF SV>>	TR FF SV>>	X5 FF SV>>	PZ FF SV>>	DW FF SV>>	TX FF SV>>	TTX FF SV>>	XX FF SV>>	PTX FF SV>>	56 FF SV>>	PC FF SV>>	PT FF SV>>	PX FF SV>>	PPZ FF SV>>	PW FF SV>>	16 FF SV>>	T6 FF SV>>	PZ6 FF SV>>	67 FF SV>>	T67 FF SV>>	X67 FF SV>>	PZ67 FF SV>>	74 FF SV>>	
Shamir Smart Office™	Power Range		-6.00 to +6.00										-6.00 to +8.00												
	Nominal Base Curve*		2, 3, 7.5, 5.25, 6.50, 8										2, 3, 5, 6.25, 8, 9.5												
	Dynamic Reduction[D]		-0.75, -1.00, -1.25, -1.50, -1.75, -2.00, -2.25										-0.75, -1.00, -1.25, -1.50, -1.75, -2.00, -2.25												
	Color	Clear											Clear												
Invisible Marking 14mm	15 FFO -175											PC FFO -175													



SHAMIR

\* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.