





Shamir Autograph Intelligence™



Minimum Fitting Height: 13mm, 15mm and 18mm

Shamir Autograph Intelligence™



Material	Invisible Marking xx	Material	Invisible Marking xx
Hard Resin 1.5	15	Polycarbonate 1.59 XTRActive®	PX
Transitions® 1.5	TR	Polycarbonate 1.59 Polarised	PPZ
XTRActive® 1.5	X5	SuperLite™ 1.6	16
Transitions® Style Colours	T	SuperLite 1.6 Transitions®	T6
Polarised 1.5	PZ	SuperLite 1.6 XTRActive®	X6
Drivewear® 1.5	DW	SuperLite 1.6 Polarised	PZ6
Vantage™ 1.5	V	Superlite 1.6 Tribrid	TriBrid
Infinite™ 1.5	INF	Superlite 1.6 Transitions® Tribrid	TriBridT
Trivex DLC™ 1.53	TX	SuperLite 1.67	67
Trivex DLC 1.53 Transitions®	TTX	SuperLite 1.67 Transitions®	T67
Trivex DLC 1.53 Drivewear®	DW	SuperLite 1.67 XTRActive®	X67
Trivex DLC 1.53 Polarised	PTX	SuperLite 1.67 Polarised	PZ67
Trivex DLC 1.53 Vantage	VX	SuperLite 1.74	74
Polycarbonate 1.59	PC	SuperLite 1.74 Transitions®	T74
Polycarbonate 1.59 Transitions®	PT		