

Shamir FirstPAL™

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Add	Color	Fitting Height	Invisible Marking	
Hard Resin	1.498	-12.25 to +6.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8	0.75 to 1.50	Clear	15mm	15 S FP	S 15
Transitions®	1.5	-12.25 to +6.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8	0.75 to 1.50	Gray & Brown	15mm	TR S FP	S 15
Transitions® XTRActive™	1.5	-12.25 to +6.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8	0.75 to 1.50	Gray	15mm	X5 S FP	S 15
Polarized	1.498	-12.25 to +6.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8	0.75 to 1.50	Gray & Brown	15mm	PZ S FP	S 15
Polycarbonate	1.586	-14.75 to +6.00	0.5, 2, 3, 4, 5.25, 6.25, 7.5	0.75 to 1.50	Clear	15mm	PC S FP	S 15
Polycarbonate Transitions®	1.586	-14.75 to +6.00	0.5, 2, 3, 4, 5.25, 6.25, 7.5	0.75 to 1.50	Gray & Brown	15mm	PT S FP	S 15
Polycarbonate Transitions® XTRActive™	1.586	-14.75 to +6.00	0.5, 2, 3, 4, 5.25, 6.25, 7.5	0.75 to 1.50	Gray	15mm	PX S FP	S 15
Polycarbonate Polarized	1.586	-14.75 to +6.00	0.5, 2, 3, 4, 5.25, 6.25, 7.5	0.75 to 1.50	Gray & Brown	15mm	PPZ S FP	S 15
Superlite™ 1.6	1.6	-15.00 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 1.50	Clear	15mm	16 S FP	S 15
Superlite 1.6 Transitions®	1.6	-15.00 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 1.50	Gray & Brown	15mm	T6 S FP	S 15
Superlite 1.6 Transitions® XTRActive™	1.6	-15.00 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 1.50	Gray	15mm	X6 S FP	S 15
Superlite 1.6 Polarized	1.6	-15.00 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 1.50	Gray & Brown	15mm	PZ6 S FP	S 15

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



SHAMIR

ReCreating Perfect Vision®

Shamir FirstPAL™

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Add	Color	Fitting Height	Invisible Marking
Superlite 1.67	1.67	-16.75 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 1.50	Clear	15mm	67 S FP 15
Superlite 1.67 Transitions®	1.67	-16.75 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 1.50	Gray & Brown	15mm	T67 S FP 15
Superlite 1.67 Transitions® XTRActive™	1.67	-16.75 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 1.50	Gray	15mm	X67 S FP 15
Superlite 1.67 Polarized	1.67	-16.75 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 1.50	Gray & Brown	15mm	PZ67 S FP 15

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



ReCreating Perfect Vision®