

Shamir Smart Attitude™

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Color	Invisible Marking
Hard Resin	1.498	-6.00 to +6.00	6, 7, 8	Clear	15 ☉ FF SA
Transitions®	1.5	-6.00 to +6.00	6, 7, 8	Gray & Brown	TR ☉ FF SA
Transitions® XTRActive™	1.5	-6.00 to +6.00	6, 7, 8	Gray	X5 ☉ FF SA
Polarized	1.498	-6.00 to +6.00	6, 7, 8	Gray & Brown	PZ ☉ FF SA
Drivewear®	1.5	-6.00 to +6.00	6, 7, 8	Drivewear®	DW ☉ FF SA
DLC™ (Trivex)	1.53	-6.00 to +6.00	6, 7, 8	Clear	TX ☉ FF SA
DLC™ (Trivex) Transitions®	1.53	-6.00 to +6.00	6, 7, 8	Gray & Brown	TTX ☉ FF SA
DLC™ (Trivex) Polarized	1.53	-6.00 to +6.00	6, 7, 8	Gray & Brown	PTX ☉ FF SA
1.56	1.56	-6.00 to +6.00	6, 7, 8	Clear	56 ☉ FF SA
Polycarbonate	1.586	-6.50 to +6.00	4, 6, 7, 8	Clear	PC ☉ FF SA
Polycarbonate Transitions®	1.67	-6.50 to +6.00	4, 6, 7, 8	Gray & Brown	PT ☉ FF SA

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



ReCreating Perfect Vision®

Shamir Smart Attitude™

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Color	Invisible Marking
Polycarbonate Transitions® XTRActive™	1.586	-6.50 to +6.00	4, 6, 7, 8	Gray	PX ☉ FF SA
Polycarbonate Polarized	1.586	-6.50 to +6.00	4, 6, 7, 8	Gray & Brown	PPZ ☉ FF SA
Polycarbonate Drivewear®	1.586	-6.50 to +6.00	4, 6, 7, 8	Drivewear®	PD ☉ FF SA
Superlite™ 1.6	1.6	-6.50 to +6.00	4, 6, 7, 8, 9, 10	Clear	16 ☉ FF SA
Superlite 1.6 Transitions®	1.6	-6.50 to +6.00	4, 6, 7, 8, 9, 10	Gray & Brown	T6 ☉ FF SA
Superlite 1.6 Transitions® XTRActive™	1.6	-6.50 to +6.00	4, 6, 7, 8, 9, 10	Gray	X6 ☉ FF SA
Superlite 1.6 Polarized	1.6	-6.50 to +6.00	4, 6, 7, 8, 9, 10	Gray & Brown	PZ6 ☉ FF SA
Superlite 1.67	1.67	-6.00 to +6.00	6, 7, 8	Clear	67 ☉ FF SA
Superlite 1.67 Transitions®	1.67	-6.00 to +6.00	6, 7, 8	Gray & Brown	T67 ☉ FF SA
Superlite 1.67 Transitions® XTRActive™	1.67	-6.00 to +6.00	6, 7, 8	Gray	X67 ☉ FF SA
Superlite 1.67 Polarized	1.67	-6.00 to +6.00	6, 7, 8	Gray & Brown	PZ67 ☉ FF SA

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



ReCreating Perfect Vision®