

Shamir Urban™

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Add	Color	Fitting Height	Invisible Marking	
Hard Resin	1.498	-12.25 to +6.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8	0.75 to 3.50	Clear	16mm and 19mm	15 URX	250
Transitions®	1.5	-12.25 to +6.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8	0.75 to 3.50	Gray & Brown	16mm and 19mm	TR URX	250
Transitions® XTRActive®	1.5	-12.25 to +6.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8	0.75 to 3.50	Gray	16mm and 19mm	X5 URX	250
Polarized	1.498	-12.25 to +6.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8	0.75 to 3.50	Gray & Brown	16mm and 19mm	PZ URX	250
Drivewear®	1.5	-12.25 to +6.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8	0.75 to 3.50	Drivewear®	16mm and 19mm	DW URX	250
DLC™ (Trivex)	1.53	-13.25 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Clear	16mm and 19mm	TX URX	250
DLC™ (Trivex) Transitions®	1.53	-13.25 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Gray & Brown	16mm and 19mm	TTX URX	250
DLC™ (Trivex) Transitions® XTRActive™	1.53	-13.25 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Gray	16mm and 19mm	TXX URX	250
DLC™ (Trivex) Polarized	1.53	-13.25 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Gray & Brown	16mm and 19mm	PTX URX	250
1.56	1.56	-14.00 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Clear	16mm and 19mm	56 URX	250
Polycarbonate	1.586	-14.75 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Clear	16mm and 19mm	PC URX	250
Polycarbonate Transitions®	1.586	-14.75 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Gray & Brown	16mm and 19mm	PT URX	250

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



ReCreating Perfect Vision®

Shamir Urban™

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Add	Color	Fitting Height	Invisible Marking
Polycarbonate Transitions® XTRActive™	1.586	-14.75 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Gray	16mm and 19mm	PX URX 250
Polycarbonate Polarized	1.586	-14.75 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Gray & Brown	16mm and 19mm	PPZ URX 250
Polycarbonate Drivewear®	1.586	-14.75 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Drivewear®	16mm and 19mm	PW URX 250
Superlite™ 1.6	1.6	-15.00 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Clear	16mm and 19mm	16 URX 250
Superlite™ 1.6 Transitions®	1.6	-15.00 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Gray & Brown	16mm and 19mm	T6 URX 250
Superlite™ 1.6 Transitions® XTRActive™	1.6	-15.00 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Gray	16mm and 19mm	X6 URX 250
Superlite™ 1.6 Polarized	1.6	-15.00 to +6.00	0.5, 2, 3, 4, 5, 6, 7	0.75 to 3.50	Gray & Brown	16mm and 19mm	PZ6 URX 250
Superlite™ 1.67	1.67	-16.75 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Clear	16mm and 19mm	67 URX 250
Superlite™ 1.67 Transitions®	1.67	-16.75 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Gray & Brown	16mm and 19mm	T67 URX 250
Superlite™ 1.67 Transitions® XTRActive™	1.67	-16.75 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Gray	16mm and 19mm	X67 URX 250
Superlite™ 1.67 Polarized	1.67	-16.75 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Gray & Brown	16mm and 19mm	PZ67 URX 250
Superlite™ 1.74	1.74	-18.75 to +15.50	0.5, 1.5, 2.5, 3.5, 4.75, 6, 7.25, 8.25, 9, 9.75, 10.5, 11.25, 12	0.75 to 3.50	Clear	16mm and 19mm	74 URX 250
Superlite™ 1.74 Transitions®	1.74	-18.75 to +15.50	0.5, 1.5, 2.5, 3.5, 4.75, 6, 7.25, 8.25, 9, 9.75, 10.5, 11.25, 12	0.75 to 3.50	Gray & Brown	16mm and 19mm	X74 URX 250

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



ReCreating Perfect Vision®